Compare Sterilizer Lifetime Cost*



Some sterilizers may seem affordable up front but demand a lifetime investment of thousands more for upkeep. While every sterilizer requires maintenance, more frequent or extensive requirements lead to higher lifetime costs. See how these sterilizers compare.



Midmark M11®

Dynamic Air Removal: Steam Flush Pressure Pulse



Tuttnauer T-Edge® 11

Dynamic Air Removal: Prevacuum (as optional upgrade)



W&H Lexa®

Dynamic Air Removal: Steam Flush Pressure Pulse with Vacuum-Assisted Drying



SciCan StatClave® G4

Dynamic Air Removal: Prevacuum

TOTAL ESTIMATED COSTS

AVERAGE SALE PRICE + LIFETIME STANDARD MAINTENANCE COSTS					
\$10,044	\$36,561	\$10,045	\$25,490		
CYCLES PER LIFETIME					
25,000 CYCLES	10,000 CYCLES	8,000 CYCLES	10,000 CYCLES		
LIFETIME COST PER CYCLE					
\$0.40	\$3.66	\$1.26	\$2.55		
TOTAL COST FOR 25,000 CYCLES					
\$10,044	\$91,403	\$31,391	\$63,725		

MAINTENANCE REQUIREMENTS

SERVICE BY TECHNICIAN RECOMMENDED					
AFTER 25,000 CYCLES OR 10 YEARS	EVERY 2 MONTHS	FREQUENCY NOT STATED; SIMILAR MODELS STATE EVERY 4,000 CYCLES OR 5 YEARS	EVERY 1,000 CYCLES OR 1 YEAR		
ANNUAL PREVENTIVE MAINTENANCE KIT					
NO CDC DAILY VACUUM TEST (BOWIE DICK) REQUIRED					
	8				
NO HEPA FILTER REQUIRED					
			×		

Helpful Tip: With reminders, step-by-step tutorials and automated history storage for routine care, Midmark helps make maintenance even easier. **Save now with the cost-effective operation and simple maintenance of Midmark sterilizers.**

*Chart data is gathered from manufacturers' statements and published materials. A lifetime is considered 10 years. Lifetime standard maintenance costs are calculated based on average prices of consumables and service fees for the required and recommended maintenance for each sterilizer over 10 years. Each total cost for 25,000 cycles is calculated based on average sterilizer selling price plus standard maintenance costs up to 25,000 cycles.